

Jazzy[®]

1450

The Jazzy[®] 1450 is a front-wheel drive bariatric power base which has a weight capacity of 272 kg (600 lbs.) and durable construction to make it an exceptionally rugged power chair. Its unique design delivers high-performance operation indoors and out.



Pride[®]
Mobility Products
Australia Pty. Ltd.

Specifications

Weight capacity	272 kg (600 lbs.)
Maximum speed	Up to 4.5 mph (7.2 km/h)
Ground clearance	3.81 cm (1.5") at anti-tip wheels
Turning radius	67.3 cm (26.5") with front riggings
Overall length	97.16 cm (38.25") without foot riggings 114.94 cm (45.25") with foot riggings
Overall width	74.3 cm (29.25")
Drive wheels	35.56 cm (14")
Front wheels	10.16 cm (4") solid anti-tips
Rear wheels	22.9 cm (9") solid castors
Suspension	Swing Arm Suspension
Drivetrain	Two-motor, Front-Wheel Drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	90 amp, PG VR2
Specialty controls	N/A
Battery charger	Off-board, 8 amp
Battery requirements	(2) 12 volt, deep cycle
Battery size	70 AH (Group 24)
Battery weight	23.6 kg (52 lbs.) each
Base weight	88.5 kg (195 lbs.)
Standard seat weight	23.4 kg (51.5 lbs.) (high-back)
Max. seat dimensions	81.3 cm x 71.1 cm (32" x 28")


The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

Features


- 272 kg (600 lbs.) weight capacity
- Maximum speed up to 4.5 mph (7.2 km/h)
- High-performance, high-efficiency motor package
- Seat dimensions up to 81.3 cm x 71.1 cm (32" x 28")

Options

- Synergy® Seat
- High-back seat with solid seat pan
- Cup holder
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder



Candy
Apple Red



Viper Blue



Pride[®]
Mobility Products
Australia Pty. Ltd.
— LIVE YOUR BEST —

20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

(Tel) (03) 8770 9600 • (Fax) (03) 9703 2720 • www.pridemobility.com.au



Quality
Endorsed
Company
9001 QEC20029
SAI Global