

SilverBoom 300



The SilverBoom 300 easily fits into most vehicles and is a perfect complement to a Jazzy® Power Chair or Pride® Scooter. Engineered to provide smooth, efficient operation, maximum durability, and sporty style, the SilverBoom 300 makes it easy and affordable to take your mobility product with you wherever you go.

Features

- Lifting capacity 136 kg (300 lbs.)
- Low headroom required - 76 cm (30")
- Boom arm end cap to protect vehicle headliners
- End swivel hook to prevent strap twisting
- Stop warning on strap to identify lifting limit
- Lift is easily disassembled for stowing
- Integrated boom arm securement when not stowed
- Remote hand control operation
- Telescoping base legs and extending boom arm length for maximum installation adjustability
- Integrated dual protective devices to prevent lift and vehicle electrical system damage

Specifications

Base weight:	9 kg (20.5 lbs) with manual base ¹
Maximum lifting capacity:	136 kg (300 lbs.)
Maximum travel distance from the ground:	203 cm (80")
Motor:	Sealed 12 volt DC
Vehicle platform:	Full-size van, Mini van, Sedan, SUV and Crossover
Warranty:	2-year transferable



Lift ready to load the product.



Loading the product.



Product loaded into vehicle.


Pride SilverStar® Lifts & Ramps

A Division of Pride® Mobility Products Australia


21 Healey Rd. • Dandenong, Victoria 3175 • www.silverstarmobility.com • www.pridemobility.com
(Tel) 03 9706 4611 • (Fax) 03 9706 4622

SilverBoom 300


Arm Height




Post Height




Base Footprint



Arm Reach



Strap Length



Can easily lift the following current Power Chairs ²

- Go-Chair™
- Jazzy® 600
- Jazzy® 1113 ATS
- Jazzy® Select™
- Jazzy® Select™ 6
- Jazzy® Select™ 14
- Jazzy® Select™ 14 XL
- Jet® 3 Ultra

Can easily lift the following current Quantum® Power Chairs ²

- Quantum® 600
- Quantum® 610
- Quantum® 1121
- Quantum® 1420

FOOTNOTES:

- ¹ Includes post weight
- ² Seating system dependent

Can easily lift the following current Pride® Scooters

- Go-Go Elite Traveller®
- Celebrity® X
- Celebrity® XL
- Hurricane®
- Legend®
- Legend® XL
- Maxima
- Revo®
- Pathrider 10

Pride SilverStar® Lifts & Ramps

A Division of Pride® Mobility Products Australia

21 Healey Rd. • Dandenong, Victoria 3175 • www.silverstarmobility.com • www.pridemobility.com

(Tel) 03 9706 4611 • (Fax) 03 9706 4622

ABN: 66 088 609 661 ACN: 088 609 661



Pride®
Mobility Products
Australia Pty. Ltd.