

QUANTUM
Jazzy

Jazzy® 1420

The Jazzy® 1420 offers an incredible balance of performance, styling and rehab capabilities. Designed especially for people with rehab needs, the 1420 offers an impressive list of standard features.

The 1420 accepts a broad range of seating options. The fully programmable digital controller uses plug-and-play technology to interface with a wide variety of specialized driver controls. In addition, the 1420's 41 cm (16") knobby drive wheels, adjustable front anti-tips, Active-Trac Suspension and patented two-motor mid-wheel drive deliver superior outdoor performance without compromising indoor manoevrability. Equipped with a standard high-torque motor package and available with all Jazzy options and accessories, the Jazzy 1420 is the ideal Power Chair for all your rehab needs.



Pride
Mobility Products
Australia Pty. Ltd.

Jazzy® 1420 shown in Candy Apple Red
with Euro Comfort seat, Pilot
Plus controller and height
adjustable footboard

Specifications

Drive Wheels:	41 cm (16") pneumatic; mid-wheel
Rear Castors:	23 cm (9") pneumatic; articulating
Anti-tip:	20 cm (8") solid; adjustable front
Suspension:	Independent rear and Active-Trac Suspension
Maximum Speed:	Adjustable up to 9 kmph (6mph), variable
Brakes:	"Intelligent Braking" electronic regenerative, disc park brake
Ground Clearance:	11 cm (4.25")
Turning Radius:	62 cm (24.5")
Overall Size:	Length: 114 cm (45") Width: 67 cm (26.5")
Standard Seating:	Euro Comfort
Seating Options:	Euro Sling Synergy Seat Deluxe High Back Versa Tilt
Drive Train:	Two motor mid-wheel drive
Battery:	Two 70 AH Group 24
Per-charge Range:	40 km (25 miles) varies with passenger weight and conditions
Battery Charger:	On-board
Electronics:	Penny and Giles 70 amp Pilot Plus controller
Weight Capacity:	170kg (375 lbs)
Standard Colours:	Candy Apple Red, Viper Blue and Forest Green
Optional Colours:	Purple, Black and Champagne
Warranty:	2-year limited warranty

Jazzy® 1420 shown in Viper Blue with Euro Sling Seat, Pilot Plus controller and swing-away footrests



Features

- Independent rear and Active-Trac suspension for outstanding outdoor performance
- Adjustable front anti-tips
- Powerful, patented two motor mid-wheel drive design provides *outstanding* manoeuvrability
- Car transportable (Not available with Versa Tilt Seat Option)
- Low centre of gravity design delivers a very smooth and stable ride
- 41 cm (16") knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 9 kmph (6 mph)
- Up to a weight capacity of 170kg (375 lbs)
- Ergonomically enhanced seating

Options

- Cup Holder
- Rear Basket
- Oxygen Tank Holder
- Cane/Crutch Holder
- Flat-free Inserts
- Swing-away Footrests
- Elevating Leg Rests
- Heel Loops
- Adjustable-Height Armrests
- Multiple Seating Options
- Specialty Controls

Jazzy® 1420 shown in Forest Green with Synergy Seat, Pilot Plus controller and swing-away footrests



Pride
Mobility Products
Australia Pty. Ltd.

21 Healey Rd. Dandenong 3175
Victoria, Australia
phone: (03) 9706-4611
fax: (03) 9706-4622
ACN: 088 609 661 ABN: 66 088 609 661

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.