



<43.5 km
Range per Charge



159 kg
Weight Capacity



158 cm
Turning Radius



More Room to Roam

Experience the ultimate in versatility and adaptability with the Celebrity MFP. Its full suspension and pillow top seating ensure a comfortable ride over various terrains. With features like a water-resistant digital display screen, rear bumper guard, a delta steering wheel, and supplemental hand brake, travelling is effortless, no matter your destination. The long range per charge and full LED lighting package make this mobility scooter your ideal travel companion, whether at the park or exploring urban cityscapes. Nothing can hold you back with the Celebrity MFP.





Rear Bumper





Digital Display



Supplemental Hand Brake

RANGE PER CHARGE (MAX) ^{1,3}	Up to 43.5 km
MAXIMUM WEIGHT CAPACITY	159 kg
TURNING RADIUS ²	158 cm
MAXIMUM SPEED ¹	10 km/h
LENGTH & WIDTH ^{2,5}	L: 145 cm W: 70.7 cm
SEAT TYPE DEPTH & WIDTH (STANDARD)	Pillow top deluxe H-B w/ sliders D: 44.5 cm (17.5") W: 49.5 cm (19.5")
SEAT-TO-DECK HEIGHT ²	45.2 cm
SEAT-TO-GROUND HEIGHT ²	66.7 cm - 70.5 cm
GROUND CLEARANCE ^{2,4}	8.8 cm at motor

OVERALL WEIGHT (NOT INC BATTERY) ^{2,6}	104.3 kg
WEIGHT OF HEAVIEST PIECE ²	94 kg (main frame does not disassemble)
SUSPENSION	Front and rear
3W/4W OPTION	4W
FRONT TYRES ²	13"
REAR TYRES ²	13"
BATTERY REQUIREMENTS ⁸	(2) 70AH
BATTERY CHARGER	Off-board, 8A
BATTERY WEIGHT ²	22 kg ⁷
COLOUR OPTIONS	 Red  Blue

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
4. Does not account for anti-tip clearance.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Durable solid plastic.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.