

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



PRODUCT CATALOGUE



Go Go Endurance Li

Leading the Way with Innovation & Quality

Dedicated to helping those with mobility needs live the highest quality of life, Pride Mobility Products® New Zealand is the New Zealand subsidiary of Pride Mobility Products Corporation, the world's leading innovator and manufacturer of personal mobility products.

At Pride, consumer needs and wishes are the driving force. We're dedicated to not just meeting medical and clinical needs, but also quality-of-life needs. It's this passion for enhancing the user experience that's propelled our growth over the years.

The Pride family of brands includes the iconic Jazzy® brand of power chairs, Go Go® Travel Mobility, Pride Scooters, and Pride Power Lift Recliners. It also encompasses Quantum Rehab®, a leader in complex rehab power chairs with iLevel® power seat elevation technology, Q6® Series of power bases, TRU-Balance® 3 and 4 Power Positioning System, and Q-Logic 3 Advanced Drive Control System, as well as Stealth Products®, the foremost manufacturer of positioning components.

With over two decades of serving New Zealand, the Pride family of brands continues creating user-inspired products with the industry's most innovative features, distinctive styles, and performance-driven touches to help consumers lead the highest quality of life possible. As has been our guiding principle, Pride's mission is to help all "Live Your Best!"

MOBILITY SCOOTERS



Quest

QUEST



Perfect for those with an active lifestyle. The compact design simply folds down in just three steps, allowing for convenient storage in a cupboard or vehicle boot. The Quest includes a lithium ion battery, which is ideal for air travel.


 Range per Charge: <19.3 km

 Weight Capacity: 113.3 kg

 Turning Radius: 82 cm

GO GO[®] LX WITH CTS SUSPENSION



 The Go Go[®] LX with Comfort-Trac Suspension (CTS) offers revolutionary technology for the smoothest, most comfortable ride. Enjoy travelling from your backyard to the park with ease.

 Range per Charge: <16.9 km

 Weight Capacity: 136 kg

 Turning Radius: 116.2 cm

GO GO[®] ENDURANCE Li

NEW



A convenient travel companion that provides an effortless mode of transportation for people who are always on the go. Easily disassembles into 5 portable pieces and includes a Lithium ion battery, allowing you to reach more destinations on a single charge.



Range per Charge: <23 km

Weight Capacity: 147 kg

Turning Radius: 122.5 cm

REVO[®] 2.0



Featuring Pride's Comfort-Trac Suspension (CTS) and feather touch disassembly, The Revo[®] 2.0 is easy to transport and ideal for those wanting a scooter with a higher weight capacity and stability for indoor and outdoor use.



Range per Charge: <21 km

Weight Capacity: 170 kg

Turning Radius: 132.1 cm

PATHRIDER[®] ES 10



With great outdoor stability, ease of use and maximum value, The Pathrider[®] ES 10 features thumb control throttle levers, a bright headlight and adjustable tiller. It also disassembles quickly offering hassle free transport.



Range per Charge: <17.9 km

Weight Capacity: 158.7 kg

Turning Radius: 132 cm



Revo 2.0





Go Go LX with CTS Suspension


PATHRIDER® 10 DELUXE

The Pathrider® 10 Deluxe offers superior comfort and stability for those seeking a safe and enjoyable ride. With features including Pride's CTS suspension, a fully adjustable high back seat and LED lighting package.



 Range per Charge: <25 km


 Weight Capacity: 158.7 kg


 Turning Radius: 133.3 cm


VICTORY® 10 LX WITH CTS SUSPENSION

Rough roads smooth out in front of you with the Victory® 10 LX, featuring Pride's front and rear Comfort-Trac suspension (CTS). Along with an LED lighting package and a deluxe high-back seat, which ensures a comfortable ride every time.



 Range per Charge: <26 km

 Weight Capacity: 181.4 kg

 Turning Radius: 152 cm



Pathrider 140 XL

CELEBRITY[®] XL DELUXE



The Celebrity XL Deluxe scooter combines a bold outdoor design with luxurious features. With an angle adjustable independent console, high-back reclining captain's seat and full directional light package, it's the ultimate outdoor scooter in its class.

 Range per Charge: <40.2 km

 Weight Capacity: 158.7 kg

 Turning Radius: 168.2 cm




Go Go Endurance Li


PATHRIDER[®] 130 XL



The Pathrider[®] 130 XL incorporates a blend of power, style and performance to handle rugged outdoor terrain. The powerful drivetrain, 24 Volt, 4-pole motor for increased power, front & rear suspension and high back reclining seat will accelerate your active lifestyle.

 Range per Charge: <25 km


 Weight Capacity: 158.7 kg

 Turning Radius: 170 cm

PATHRIDER[®] 140 XL



Discover the perfect blend of power and precision with the Pathrider[®] 140 XL. Equipped with luxurious features, such as front and rear suspension, a high-back pillow-top seat, and an LED lighting package, it offers the perfect outdoor mobility solution.

 Range per Charge: <35 km

 Weight Capacity: 181.4 kg

 Turning Radius: 208 cm

OUTBACK®

Rugged and built to last, the Outback® is your perfect adventure companion! Equipped with front and rear suspension, LED lighting, a user-friendly console and dual hydraulic brakes, the Outback® is designed for embracing the great outdoors!



 Range per Charge: <36 km


 Weight Capacity: 158.7 kg

 Turning Radius: 243 cm



Outback

FOR MORE INFORMATION
VISIT OUR WEBSITE




www.pridemobility.co.nz

SPORT RIDER

The Sport Rider is an extraordinary three-wheel motorcycle inspired outdoor scooter. This scooter boasts impressive features including a full digital display for speed and distance tracking, a robust suspension system, an adjustable high-back seat, and a powerful lighting kit.

 Range per Charge: <35 km

 Weight Capacity: 181.4 kg

 Turning Radius: 158 cm



END-TO-END | Scooter Specifications



	QUEST®	GO GO® LX W/ CTS SUSPENSION	GO GO ENDURANCE LI®	REVO® 2.0	PATHRIDER® ES10	PATHRIDER® 10 DX	VICTORY® 10 LX W/ CTS SUSPENSION	CELEBRITY® XL DX	PATHRIDER® 130 XL	PATHRIDER® 140 XL	OUTBACK®	SPORT RIDER
RANGE PER CHARGE (MAX) ^{2, 8}	Up to 19.3 km	Up to 16.9 km ⁹	Up to 23 km	Up to 21 km	Up to 17.9 km	Up to 25 km with 40AH	Up to 26 km	Up to 40.2 km ¹⁵	Up to 25 km	Up to 35 km	Up to 35 km	Up to 35 km
MAXIMUM WEIGHT CAPACITY	113.3 kg	136 kg	147 kg	170 kg	158.7 kg	158.7 kg	181.4 kg	158.7 kg	158.7 kg	181.4 kg	158.7 kg	181.4 kg
TURNING RADIUS ³	82 cm	116.2 cm	122.5 cm	132.1 cm	132 cm	133.3 cm	152 cm	168.2 cm	170 cm	210 cm right 208 cm left	243 cm	158 cm
MAXIMUM SPEED ¹⁷	Up to 6 km/h	Up to 7 km/h	Up to 6.7 km/h	Up to 8 km/h	Up to 8.85 km/h	Up to 8.85 km/h	Up to 9 km/h	Up to 15 km/h	Up to 15 km/h	Up to 15 km/h	Up to 15 km/h	Up to 15 km/h
LENGTH & WIDTH ^{3, 11}	L: 92.7 cm W: 49.5 cm	L: 104 cm W: 52.1 cm	L: 107 cm W: 52 cm	L: 119.4 cm W: 52.7 cm	L: 115.6 cm W: 57.1 cm	L: 119.4 cm W: 56.2 cm	L: 119.4 cm W: 56.5 cm	L: 130.2 cm W: 64.1 cm	L: 136 cm W: 63.6 cm	L: 142.8 cm W: 69.2 cm	L: 156 cm W: 83 cm	L: 165 cm W: 78 cm
SEAT TYPE (STANDARD) ⁵	Compact, Foldable, Molded Plastic, 41 cm x 33 cm (16"x13")	Compact, Foldable, 43 cm x 43 cm (17" x 17")	Compact, Foldable, 43.1 cm x 43.1 cm (17" x 17") 50.8 cm x 45.7 cm (20" x 18")	Compact, Foldable, 43 cm x 43 cm (17" x 17")	Foldable Molded Plastic 46 cm x 43 cm (18" x 17")	CRS Junior Seat 44.25 cm x 43 cm (17.5" x 17")	Reclining H-B w/ Sliders & Headrest 46 cm x 46 cm (18" x 18")	Deluxe Reclining H-B w/ Sliders & Headrest 46 cm x 46 cm (18" x 18")	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Captains High-Back Adjustable Seat 51 cm x 56 cm (20" x 22")
SEAT-TO-DECK HEIGHT RANGE	31 cm Fixed	41.3 cm - 43.8 cm	41.2 cm x 45 cm	38.1 cm - 45.7 cm	37 cm - 42 cm	43.8 cm - 47.6 cm	46 cm - 49 cm	43.2 cm x 53.3 cm	42 cm - 47 cm	47 cm	50 cm - 53 cm	47.5 cm
SEAT-TO-GROUND HEIGHT RANGE	51 cm Fixed	52.7 cm - 56.5 cm	53.9 cm x 57.7 cm	53.3 cm - 60 cm	51.4 cm - 55.2 cm	60 cm - 63.8 cm	60 cm - 64.8 cm	54 cm - 64 cm	56.7 cm	66.7 cm - 70.5 cm	76 cm - 80 cm	24.8 cm
GROUND CLEARANCE ^{3, 12, 14}	5.8 cm at frame	6.4 cm at frame 3.8 cm at motor	3.6 cm at motor	8.9 cm at frame 3.2 cm at motor	8.8 cm at frame	8.25 cm at frame	5.6 cm at motor	10.8 cm at transaxle	9 cm at transaxle	13.3 cm at midframe	14 cm at midframe	12.7 cm at shroud
OVERALL WEIGHT (NOT INC. BATTERY) ^{3, 13}	23.5 kg	39.4 kg	42.4 kg	60 kg	59.4 kg	84.8 kg	81.2 kg	112 kg	99.8 kg	117.7 kg	124 kg	115 kg
WEIGHT OF HEAVIEST PIECE ³	23.5 kg (does not disassemble)	16.6 kg (front section)	18.48 kg (front section)	25.9 kg (front section)	25.8 kg (front section)	29.5 kg (front section)	35.2 kg (rear section)	37.6 kg (front section)	34 kg (rear section)	93 kg (main frame)	99 kg (main frame)	92.3 kg (main frame)
SUSPENSION	N/A	CTS Front & Rear	CTS Front & Rear	Front & Rear	N/A	Front & Rear	Front & Rear	Front & Rear	Front & Rear	Front & Rear	Front & Rear	Front & Rear
3W/4W OPTION	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	3W
FRONT TYRES	1.5" x 7" solid	2" x 7" solid	2" x 8" solid	2.5" x 9" solid	3.6" x 10.4"	3.5" x 10"	3.5" x 10"	4" x 13"	4" x 13"	4" x 13"	14.5" x 6"	2.75" x 19" ¹⁶
REAR TYRES	1.9" x 7.8" solid	2.5" x 8" solid	2.5" x 8" solid	2.5" x 9" solid	3.6" x 10.75" solid	4" x 10"	3.5" x 10"	4" x 10"	4" x 13"	4" x 13"	14.5" x 6"	2.75" x 19" ¹⁶
BATTERY REQUIREMENTS ⁶	(1) 11.5 AH (Lithium Ion)	(2) 18AH	16AH (Lithium Ion)	(2) 31AH	(2) 31AH	(2) 31AH or (2) 40AH	(2) 40AH	(2) 40AH	(2) 55AH	(2) 70AH	(2) 100AH	(2) 70AH
BATTERY CHARGER	Off-Board 2.5A (Lithium Ion)	Off-Board, 2A	Off-board, 3.5A (Lithium Ion)	Off-Board, 3.5A	Off-Board, 3A	Off-Board, 3A (32AH) Off-Board, 5A (40AH)	Off-Board, 5A	Off-Board, 5A	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A
BATTERY WEIGHT ¹⁰	2.3 kg	6 kg ⁷	5.4 kg	11 kg	10.5 kg	10.5 kg or 14 kg	14 kg	14 kg	16.8 kg	23.45 kg	27.9 kg	23.7 kg
COLOUR OPTIONS	White	Red, Blue ¹	Garnet Red, Euro Grey	Grey Street, True Blue	Red	Viper Blue	Candy Apple Red	Candy Apple Red, Silver	Candy Apple Red, Mercury Silver	Candy Apple Red, Yellow, Silver	Stone Grey	Black

¹ All shroud panels included. ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%. ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. ⁴ 12 AH battery pack weighs 8.8 kg (19.5 lbs.). ⁵ Durable solid plastic. ⁶ Deep cycle (AGM or gel cell type recommended) excluding Lithium-ion models. ⁷ 18 AH battery pack weighs 13.4 kg (29.5 lbs.). ⁸ Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. ⁹ Up to 14 miles per charge on 18 AH batteries. ¹⁰ Battery weight may vary based on manufacturer. ¹¹ Width of base (not including seat). ¹² Does not account for anti-tip clearance. ¹³ Includes standard seat weight. ¹⁴ At midframe ¹⁵ Range tested with 31-36 AH Batteries. Standard larger battery size would increase range. ¹⁶ 6.9 cm x 35.5 cm (2.75" x 14") Rim Only ¹⁷ Varies with user weight, terrain type, battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%. **NOTE** The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



Jazzy®

POWER CHAIRS



Jazzy Air 2

Litestream™ CARBON

NEW



With a strong carbon fibre frame, the Litestream® Carbon Power Wheelchair is a reliable travel companion. It folds down in just a few quick steps, making it a super lightweight effortless option to take when travelling on planes, trains and cars.

 Range per Charge: <15 km

 Weight Capacity: 136 kg

 Turning Radius: 120 cm

Jazzy® PASSPORT



Embrace the ease of lightweight travel with the Jazzy® Passport! The compact design effortlessly folds in a few simple steps, making it a breeze to transport. With convenient storage capabilities, the Jazzy® Passport ensures hassle-free traveling.

 Range per Charge: <15.4 km

 Weight Capacity: 113.3 kg

 Turning Radius: 61.5 cm



Jazzy Passport

Jazzy® SELECT 6



Experience the sporty and sophisticated design of the Jazzy® Select 6 power wheelchair. It provides powerful performance and stability with its superior suspension system. Featuring an adjustable high-back seat with manual recline.

 Range per Charge: <24.1 km

 Weight Capacity: 136 kg

 Turning Radius: 55.5 cm

GOCHAIR®



The Go Chair® NG is #1 for portability! With Pride's feather-touch disassembly, and a tight turning radius, the Go Chair® NG is ideal for indoor use and smooth surfaces, allowing you to enjoy greater independence.




 Range per Charge: <14 km

 Weight Capacity: 136 kg

 Turning Radius: 64.5 cm

DOWNLOAD OUR PQ-365 APP FOR EVERYTHING PRIDE AND QUANTUM

PQ-365



Download on the
App Store

GET IT ON
Google Play

Jazzy Air₂



Take your life to new heights with the Jazzy Air[®] 2! With a simple switch, this power wheelchair elevates 30.4 cm in just 11 seconds. Allowing you to safely drive at speeds of up to 6.4 km/h while elevated, enabling face-to-face social interactions.

 Range per Charge: <27 km

 Weight Capacity: 136 kg

 Turning Radius: 61.5 cm

Jazzy[®] ELITE HD





Experience the ultimate heavy-duty electric wheelchair with the Jazzy[®] Elite HD. This power wheelchair excels on various terrains, providing superior performance and control with its dual-motor, front-wheel drive technology, that offers a tight turning radius.

 Range per Charge: <15.7 km

 Weight Capacity: 204 kg

 Turning Radius: 61.9 cm



Stealth Products[®] Cushions feature patented Coolcore[®] technology, a chemical free, four-way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. Our cushions also use a colour guide that identifies the suggested function family by a unique patch.

FEATURES

- Chemical free
- Regulates temperature
- Manages and distributes your body heat

[CUSHIONS.STEALTHPRODUCTS.COM/](https://cushions.stealthproducts.com/)

Learn about our Stealth Products range of cushions featuring our exclusive Coolcore[®] Technology.



Scan here to access the
Stealth Cushions Catalogue



END-TO-END | Powerchair Specifications




	LITESTREAM CARBON	JAZZY® PASSPORT	GO CHAIR® NG	JAZZY® SELECT 6	JAZZY® AIR 2	JAZZY® ELITE HD
RANGE PER CHARGE (MAX) ^{9,5}	Up to 15 km	Up to 15.4 km	Up to 14km	Up to 24.1 km	Up to 27 km	Up to 15.7 km
MAXIMUM WEIGHT CAPACITY	136 kg	113.3 kg	136 kg	136 kg	136 kg	204.1 kg ⁶
TURNING RADIUS ¹	120 cm	61.5 cm	64.5 cm	55.5 cm	61.5 cm	61.9 cm
MAXIMUM SPEED ³	Up to 6 km/h	Up to 5.7 km/h	Up to 5.9 km/h	Up to 6.4 km/h	Up to 6.4 km/h	Up to 6.4 km/h
LENGTH & WIDTH ¹	L: 93 cm W: 57.7 cm	L: 101.6 cm W: 59.6 cm	L: 83.8 cm W: 57 cm	L: 90 cm W: 59.6 cm	L: 106 cm W: 61.5 cm	L: 90 cm W: 61.5 cm
WEIGHT OF HEAVIEST PIECE ⁴	16 kg	24.4 kg	16.3 kg	44.1 kg	91.6 kg ⁸	44.9 kg
GROUND CLEARANCE ¹	5.5 cm at motor	7.6 cm at motor	4.5 cm at frame	4.4 cm at motor mount	3.6 cm at motor mount	7.6 cm at footboard mount 3.8 cm at anti-tip wheel
SEAT SIZE (STANDARD)	43 cm x 48 cm (17" x 18.9")	41.9 cm x 44.4 cm (16.5" x 17.5")	43.1 cm x 45.7 cm (17" x 18")	45.7 cm x 45.7 cm (18" x 18")	40.6 cm x 45.7 cm (16" x 18") 45.7 cm x 45.7 cm (18" x 18") 50.8 cm x 45.7 cm (20" x 18")	45.7 cm x 45.7 cm (18" x 18") 50.8 cm x 50.8 cm (20" x 20") 55.8 cm x 50.8 cm (22" x 20")
SEAT WEIGHT (STANDARD)	N/A (seat is integrated with the base)	N/A (seat is integrated with the base) (sling)	13.6 kg	19.9 kg (comfort, high-back)	23 kg (high-back)	23.5 kg (high-back)
HIGH-BACK SEAT-TO-FLOOR RANGE	N/A	N/A	N/A	53.9 cm - 59 cm	53.9 cm - 84 cm	50.8 cm - 55.8 cm
SPECIALTY SEATING OPTIONS	N/A	N/A	N/A	N/A	N/A	N/A
ELECTRONICS (STANDARD)	35A, Micon	50A, PG VSI	40A, PG GC3	40A, PG GC3	40A, VR2	60A, Dynamic Linx
SPECIALTY CONTROLS	N/A	N/A	N/A	N/A	N/A	N/A
DRIVE TRAIN	2 x 24 V/250W with brushless motor (aluminium alloy material)	2 x 180W brushed motors rear-wheel-drive	Two-motor, rear-wheel drive	Two-motor, in-line, mid-wheel drive	Two-motor, in-line, mid-wheel drive	Two-motor, front-wheel drive
BRAKING SYSTEM	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical
SUSPENSION TYPE	Front	Limited	Limited	Active-Trac® (ATX)	Active-Trac® with Mid-Wheel 6 technology	Limited
FRONT WHEELS	7" solid castors	8" solid castors	6" solid castors	6" solid castors	6" solid castors	3" solid anti-tips
DRIVE WHEELS	8.5" solid castors	12" pneumatic	9" solid	10" solid	10" solid	14" solid
REAR WHEELS	N/A	N/A	3" solid anti-tips	6" solid castors	6" solid castors	6" solid castors
BATTERY REQUIREMENTS	24 V, 12 AH Lithium-ion (1)	24 V, 18 AH Lithium-ion (1)	12 V, 18 AH (2)	12 V, 35 AH (2)	12 V, 40 AH (2)	12 V, 55 AH (2)
BATTERY CHARGER	24V 2A, Lithium	2A, off-board, Lithium	2A, off-board	3.5A, off-board	5A, off-board	5A, off-board
BATTERY WEIGHT ⁷	1.8 kg	3.9 kg	14 kg battery box	10.5 kg each	14 kg each	16.8 kg each
COLOUR OPTIONS	Carbon	Pewter	Ruby Red, Pearl White	Ruby Red	Black Oynx, Pearl White, Pink Topaz, Garnet Red	Jazzy Red, Jazzy Blue

¹ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. ² Without front riggings. ³ Varies with user weight, terrain, battery condition and tire pressure. ⁴ Does not include seat or battery weights. ⁵ Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, - 5%. ⁶ Maximum weight capacity is reduced to 400 lbs. (181 kg) on 46 cm (18") deep seat configurations. ⁷ Battery weight may vary based on manufacturer. ⁸ Seat weight included (because tools are required to remove the seat). Battery weight not included. ⁹ Varies with user weight, terrain, battery type, battery condition and tire pressure. *International configuration shown **SEATING: A** For Synergy® Seating Contact your nearest Pride® stocklist for more information **B** Quantum® dealers only **NOTE** The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

PRIDE RAMPERS


SINGLE FOLD



Single Fold Ramps offer a reliable solution for transitioning curbs, thresholds, and other obstacles. They are 75 cm wide and are available in several different lengths. These high-traction ramps are portable and lightweight, making them ideal for home and travel.

 Weight Capacity: 272 kg

RUBBER THRESHOLD




Rubber Threshold Ramps provide a safe, durable means of navigating uneven household surfaces, including small steps and doorways. A hidden dovetail locking system allows for easy installation and incremental ramp height adjustment of up to 10.16cm (4"). A modular design, slip-resistant and maintenance free material makes our Rubber Threshold Ramps a practical and affordable home accessibility solution.

 Weight Capacity: 385 kg

ACCESSORIES



Sun Canopy & Flag



Rear View Mirror



Rear Pride Scooter Bag



Saddle Pride Bag



Walking Frame Holder



Oxygen Tank Holder



Rear Pod



XLR USB Charger



Mobile Phone Holder



Foldable Cup Holder



Tel: 0800 828 033

Visit pridemobility.co.nz for further details



NZ Distribution Centre: Factory 5/208 Swanson Rd. Henderson, Auckland NZ 0610

NOTE: All specifications subject to change without notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.