

## Maximise your potential

Embrace independence and reclaim your mobility. The all new Maxima mobility scooter is the perfect companion for those seeking freedom, convenience, and style. Whether you're exploring your neighbourhood, shopping or attending social gatherings, the Maxima scooter ensures you can do it all, with ease. Available in 3 and 4 wheel configuration, this ultra-heavy-duty scooter boasts a remarkable 227 kg weight capacity and includes side bumpers, a supplemental handbrake and USB charging port.











Side Bumpers

**USB Charging Point** 

3-wheel option available

RANGE PER CHARGE (MAX) <sup>1,3</sup>	Up to 31.4 km
MAXIMUM WEIGHT CAPACITY	227 kg
TURNING RADIUS <sup>2</sup>	134.6 cm (4-wheel) 124.5 cm (3-wheel)
MAXIMUM SPEED <sup>1</sup>	9.3 km/h
LENGTH & WIDTH <sup>2,5</sup>	L: 120.7 cm W: 68.6 cm
SEAT TYPE (STANDARD)	Adjustable High Back with Black Vinyl L: 50.8 cm (20") W: 55.9 cm (22")
SEAT-TO-DECK HEIGHT <sup>2</sup>	43.5 cm - 48.6 cm
SEAT-TO-GROUND HEIGHT <sup>2</sup>	56.1 cm - 61.2 cm
GROUND CLEARANCE <sup>2,4</sup>	5.7 cm at motor

OVERALL WEIGHT (NOT INC BATTERY) <sup>2,6</sup>	83.2 kg (3-wheel) 89.3 kg (4-wheel)
WEIGHT OF HEAVIEST PIECE <sup>2</sup>	N/A (does not disassemble)
SUSPENSION	N/A
3W/4W OPTION	3W and 4W
FRONT TYRES <sup>2</sup>	10" x 3.6" Solid
REAR TYRES <sup>2</sup>	10.75" x 3.6" Solid
BATTERY REQUIREMENTS <sup>8</sup>	(2) 50AH
BATTERY CHARGER	Off-board, 5A
BATTERY WEIGHT <sup>2</sup>	16.6 kg
COLOUR OPTIONS	Garnet Red Ocean Blue

- Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
- 2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
- 4. Does not account for anti-tip clearance.
- 5. Width of base (not including seat).
- 6. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.





