

PRIDE®
M O B I L I T Y
L I V E Y O U R B E S T ®



PRODUCT CATALOGUE



Go Go Endurance Li

Leading the Way with Innovation & Quality

Dedicated to helping those with mobility needs live the highest quality of life, Pride Mobility Products® Australia is the Australian subsidiary of Pride Mobility Products Corporation, the world's leading innovator and manufacturer of personal mobility products.

At Pride, consumer needs and wishes are the driving force. We're dedicated to not just meeting medical and clinical needs, but also quality-of-life needs. It's this passion for enhancing the user experience that's propelled our growth over the years.

The Pride family of brands includes the iconic Jazzy® brand of power chairs, Go Go® Travel Mobility, Pride Scooters, and Pride Power Lift Recliners. It also encompasses Quantum Rehab®, a leader in complex rehab power chairs with iLevel® power seat elevation technology, Q6® Series of power bases, TRU-Balance® 3 and 4 Power Positioning System, and Q-Logic 3 Advanced Drive Control System, as well as Stealth Products®, the foremost manufacturer of positioning components.

With over two decades of serving Australia, the Pride family of brands continues creating user-inspired products with the industry's most innovative features, distinctive styles, and performance-driven touches to help consumers lead the highest quality of life possible. As has been our guiding principle, Pride's mission is to help all "Live Your Best!"

MOBILITY SCOOTERS FOR TRAVEL



Quest

QUEST



Perfect for those with an active lifestyle. The compact design simply folds down in just three steps, allowing for convenient storage in a cupboard or vehicle boot. The Quest includes a lithium ion battery, which is ideal for air travel.

 Range per Charge: <19.3 km

 Weight Capacity: 113.3 kg

 Turning Radius: 82 cm

GO-GO[®] ULTRAX



Enjoy worry-free travel with the Go-Go[®] Ultra X, an ultra-compact scooter that is simple to disassemble. Available in three or four-wheel configuration.

 Range per Charge: <11.5 km

 Weight Capacity: 118 kg

 Turning Radius: 100.3 cm

GO GO[®]

ELITE TRAVELLER[®] PLUS



Perfect for travelling, PLUS so much more! The Go Go Elite Traveller[®] Plus keeps you moving through all of life's adventures, while offering easy disassembly, an ergonomic tiller and non-scuffing tyres.



Range per Charge: <20.9 km

Weight Capacity: 136 kg

Turning Radius: 116.2 cm

GO GO[®]

LX WITH CTS SUSPENSION



The Go Go[®] LX with Comfort-Trac Suspension (CTS) offers revolutionary technology for the smoothest, most comfortable ride. Enjoy travelling from your backyard to the park with ease.



Range per Charge: <16.9 km

Weight Capacity: 136 kg

Turning Radius: 116.2 cm

GO GO[®]

ENDURANCE Li



A convenient travel companion that provides an effortless mode of transportation for people who are always on the go. Easily disassembles into 5 portable pieces and includes a Lithium ion battery, allowing you to reach more destinations on a single charge.



Range per Charge: <23.1 km

Weight Capacity: 147 kg

Turning Radius: 122.5 cm

ZERO

TURN 8



Featuring Prides *iTurn Technology*[®], the Zero Turn 8 is a 4-wheel scooter with an extremely tight turning radius and features Comfort-Trac Suspension (CTS) to ensure a smooth road ahead. This scooter will have you navigating tight corners and small spaces effortlessly.



Range per Charge: <16 km

Weight Capacity: 136 kg

Turning Radius: 97.2 cm

NEW

END-TO-END | Travel Scooter Specifications



	QUEST®	GO GO® ULTRA X	GO GO® ELITE TRAVELLER PLUS	GO GO® LX W/ CTS SUSPENSION	GO GO ENDURANCE LI®	ZERO TURN 8
RANGE PER CHARGE (MAX) ^{2, 8}	Up to 19.3 km	3W: Up to 11 km 4W: Up to 11.5 km	Up to 20.9 km ⁹	Up to 16.9 km ⁹	Up to 23 km	Up to 16 km
MAXIMUM WEIGHT CAPACITY	113.3 kg	118 kg	136 kg	136 kg	147 kg	136 kg
TURNING RADIUS ³	82 cm	3W: 83.8 cm 4W: 100.3 cm	116.2 cm	116.2 cm	122.5 cm	97.2 cm
MAXIMUM SPEED ¹⁷	Up to 6 km/h	Up to 6.4 km/h	Up to 7.2 km/h	Up to 7 km/h	Up to 6.7 km/h	Up to 6.4 km/h
LENGTH & WIDTH ^{3, 11}	L: 92.7 cm W: 49.5 cm	3W: L: 94 cm 3W W: 49.5 cm 4W: L: 100.3 cm 4W W: 49.5 cm	L: 102.3 cm W: 54 cm	L: 104 cm W: 52.1 cm	L: 107 cm W: 52 cm	L: 116 cm W: 50.8 cm
SEAT TYPE (STANDARD) ⁵	Compact, Foldable, Molded Plastic, 41 cm x 33 cm (16"x13")	Compact, Foldable, 43 cm x 43 cm (17" x 17")	Compact, Foldable, 43 cm x 43 cm (17" x 17")	Compact, Foldable, 43 cm x 43 cm (17" x 17")	Compact, Foldable 43.1 cm x 43.1 cm (17" x 17") 50.8 cm x 45.7 cm (20" x 18")	Compact, Foldable, 43 cm x 43 cm (17" x 17")
SEAT-TO-DECK HEIGHT RANGE	31 cm Fixed	40 cm - 43.8 cm	41.3 cm - 45 cm	41.3 cm - 43.8 cm	41.2 cm x 45 cm	40.6 cm - 45.6 cm
SEAT-TO-GROUND HEIGHT RANGE	51 cm Fixed	55 cm - 60 cm	52 cm - 55.88 cm	52.7 cm - 56.5 cm	53.9 cm x 57.7 cm	53.3 cm - 60 cm
GROUND CLEARANCE ^{3, 12, 14}	5.8 cm at frame	6.3 cm at midframe 3.8 cm at motor	7.62 cm at frame	6.4 cm at frame 3.8 cm at motor	3.6 cm at motor	8.9 cm at frame 3.3 cm at motor
OVERALL WEIGHT (NOT INC. BATTERY) ^{3, 13}	23.5 kg	3W: 34.5 kg 4W: 36.3 kg	39.4 kg	39.4 kg	42.4 kg	53.9 kg
WEIGHT OF HEAVIEST PIECE ³	23.5 kg (does not disassemble)	3W: 13.4 kg 4W: 15.4 kg (front sections)	15.4 kg (front section)	16.6 kg (front section)	18.48 kg (front section)	21.8 kg (rear section)
SUSPENSION	N/A	N/A	N/A	CTS Front & Rear	CTS Front & Rear	CTS Front & Rear
3W/4W OPTION	4W	3W & 4W	4W	4W	4W	4W
FRONT TYRES	1.5" x 7" solid	3W: 2" x 8" solid 4W: 2" x 7" solid	2" x 8" solid	2" x 7" solid	2" x 8" solid	2" x 7" solid
REAR TYRES	1.9" x 7.8" solid	2.5" x 8" solid	3" x 9" solid	2.5" x 8" solid	2.5" x 8" solid	2.5" x 9" solid
BATTERY REQUIREMENTS ⁶	(1) 11.5 AH (Lithium Ion)	(2) 12AH	(2) 18AH	(2) 18AH	16AH (Lithium Ion)	(2) 22AH
BATTERY CHARGER	Off-Board 2.5A (Lithium Ion)	Off-Board, 2A	Off-Board, 2A	Off-Board, 2A	Off-board, 3.5A (Lithium Ion)	Off-Board, 3.5A
BATTERY WEIGHT ¹⁰	2.3 kg	4 kg ⁴	6 kg ⁸	6 kg ⁷	5.4 kg	6 kg ⁷
COLOUR OPTIONS	White	Red	Red, Blue ¹	Red, Blue ¹	Garnet Red, Euro Grey	Fire Opal, Pearl White, Sapphire Blue

1 All shroud panels included. **2** Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%. **3** Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. **4** 12 AH battery pack weighs 8.8 kg (19.5 lbs.). **5** Durable solid plastic. **6** Deep cycle (AGM or gel cell type recommended) excluding Lithium-ion models. **7** 18 AH battery pack weighs 13.4 kg (29.5 lbs.). **8** Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. **9** Up to 14 miles per charge on 18 AH batteries. **10** Battery weight may vary based on manufacturer. **11** Width of base (not including seat). **12** Does not account for

anti-tip clearance. **13** Includes standard seat weight. **14** At midframe **15** Range tested with 31-36 AH Batteries. Standard larger battery size would increase range. **16** 6.9 cm x 35.5 cm (2.75"x 14") Rim Only **17** Varies with user weight, terrain type, battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%. **NOTE** The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



MOBILITY SCOOTERS FOR EVERYDAY



Zero Turn 8

REVO® 2.0



Featuring Pride's Comfort-Trac Suspension (CTS) and feather touch disassembly, The Revo® 2.0 is easy to transport and ideal for those wanting a scooter with a higher weight capacity and stability for indoor and outdoor use.



 Range per Charge: <21 km

 Weight Capacity: 170 kg

 Turning Radius: 132.1 cm

BAJA BANDIT

NEW



The first member of the Baja® Mobility family. The Revo Bandit is bold, adventurous and innovative. With bumpers and a high ground clearance the Baja® Bandit tackles obstacles with ease and it's long-lasting lithium ion battery lets you explore for longer.



 Range per Charge: <48 km

 Weight Capacity: 181.4 kg

 Turning Radius: 132 cm

ZERO TURN 10



The Zero Turn 10 uses *iTurn Technology*[®], allowing it to turn on its own axis. Combined with Comfort-Trac Suspension (CTS) and dual motors, this scooter provides a stable and comfortable ride, indoors or outdoors.

 Range per Charge: <28.8 km

 Weight Capacity: 170 kg

 Turning Radius: 109 cm

PATHRIDER[®] ES 10



With great outdoor stability, ease of use and maximum value, The Pathrider[®] ES 10 features thumb control throttle levers, a bright headlight and adjustable tiller. It also disassembles quickly offering hassle free transport.

 Range per Charge: <17.9 km

 Weight Capacity: 158.7 kg

 Turning Radius: 132 cm

PATHRIDER[®] 10



The sleek and sporty Pathrider[®] 10 offers comfort and performance for your active lifestyle. Enjoy simple disassembly with one hand, a bright LED headlight and swivel seat to help you get through your daily activities.

 Range per Charge: <25 km

 Weight Capacity: 158.7 kg

 Turning Radius: 137 cm

PATHRIDER[®] 10 DELUXE



The Pathrider[®] 10 Deluxe offers superior comfort and stability for those seeking a safe and enjoyable ride. With features including Pride's CTS suspension, a fully adjustable high back seat and LED lighting package.

 Range per Charge: <25 km

 Weight Capacity: 158.7 kg

 Turning Radius: 133.3 cm




Revo 2.0

GOPHER[®]SV



Whether you're zipping around the city or cruising your neighbourhood, you can trust the Gopher SV to get you where you need to go. Boasting a premium lighting package, front and rear suspension and deluxe high back seat with ample leg room.

 Range per Charge: <24 km

 Weight Capacity: 160 kg

 Turning Radius: 143 cm

VICTORY[®]10 LX WITH CTS SUSPENSION



Rough roads smooth out in front of you with the Victory[®] 10 LX, featuring Pride's front and rear Comfort-Trac suspension (CTS). Along with an LED lighting package and a deluxe high-back seat, which ensures a comfortable ride every time.

 Range per Charge: <26 km

 Weight Capacity: 181.4 kg

 Turning Radius: 152 cm



CELEBRITY[®]DX



The Celebrity[®] DX combines style, performance, comfort, and includes a deluxe seat with headrest and manual recline. It's a functional mobility solution with a touch of luxury. Available in a 3 or 4-wheel configuration.

 Range per Charge: <20.6 km

 Weight Capacity: 158.7 kg

 Turning Radius: 149.8 cm

CELEBRITY[®] XL DELUXE



The Celebrity XL Deluxe scooter combines a bold outdoor design with luxurious features. With an angle adjustable independent console, high-back reclining captain's seat and full directional light package, it's the ultimate outdoor scooter in its class.



 Range per Charge: <40.2 km

 Weight Capacity: 158.7 kg

 Turning Radius: 168.2 cm

MAXIMA



Heavy duty and robust, the ALL NEW Maxima scooter maximises your potential! Experience the comfort of a new high-back, depth adjustable recline seat that swivels for easier access, side bumpers, and a USB charger standard. Available in 3 or 4-wheel configuration.



 Range per Charge: <31.4 km

 Weight Capacity: 226.7 kg

 Turning Radius: 134.6 cm

PATHRIDER[®] 130 XL



The Pathrider[®] 130 XL incorporates a blend of power, style and performance to handle rugged outdoor terrain. The powerful drivetrain, 24 Volt, 4-pole motor for increased power, front & rear suspension and high back reclining seat will accelerate your active lifestyle.



 Range per Charge: <25 km

 Weight Capacity: 158.7 kg

 Turning Radius: 170 cm

PATHRIDER[®] 140 XL



Discover the perfect blend of power and precision with the Pathrider[®] 140 XL. Equipped with luxurious features, such as front and rear suspension, a high-back pillow-top seat, and an LED lighting package, it offers the perfect outdoor mobility solution.



 Range per Charge: <35 km

 Weight Capacity: 181.4 kg

 Turning Radius: 208 cm



NEW



The Celebrity MFP offers versatility and adaptability for individuals seeking a customizable mobility scooter. It boasts a digital display screen, front and rear suspension, a steering wheel, and the ability to adjust the floor length.

 Range per Charge: <43.5 km

 Weight Capacity: 159 kg

 Turning Radius: 158 cm

OUTBACK®




Rugged and built to last, the Outback® is your perfect adventure companion! Equipped with front and rear suspension, LED lighting, a user-friendly console and dual hydraulic brakes, the Outback® is designed for embracing the great outdoors!

 Range per Charge: <35 km

 Weight Capacity: 158.7 kg

 Turning Radius: 243 cm

FOR MORE INFORMATION
VISIT OUR WEBSITE




www.pridemobility.com.au

SPORT RIDER



The Sport Rider is an extraordinary three-wheel motorcycle inspired outdoor scooter. This scooter boasts impressive features including a full digital display for speed and distance tracking, a robust suspension system, an adjustable high-back seat, and a powerful lighting kit.

 Range per Charge: <35 km

 Weight Capacity: 181.4 kg

 Turning Radius: 158 cm

END-TO-END | Scooter Specifications



	REVO® 2.0	ZERO TURN 10	PATHRIDER® ES10	PATHRIDER® 10	PATHRIDER® 10 DX	GOPHER SV	VICTORY® 10 LX W/ CTS SUSPENSION	CELEBRITY® DX	CELEBRITY® XL DX	MAXIMA
RANGE PER CHARGE (MAX) ^{2, 8}	Up to 21 km	Up to 28.8 km	Up to 17.9 km	Up to 25 km	Up to 25 km with 40AH	Up to 24 km	Up to 26 km	Up to 20.6 km	Up to 40.2 km ¹⁵	3W: Up to 31.4 km 4W: Up to 31.4 km
MAXIMUM WEIGHT CAPACITY	170 kg	170 kg	158.7 kg	158.7 kg	158.7 kg	160 kg	181.4 kg	158.7 kg	158.7 kg	226.8 kg
TURNING RADIUS ³	132.1 cm	109 cm	132 cm	137 cm	133.3 cm	143 cm	152 cm	3W: 105.4 cm 4W: 149.8 cm	168.2 cm	3W: 124.5 cm 4W: 134.6 cm
MAXIMUM SPEED ¹⁷	Up to 8 km/h	Up to 10 km/h	Up to 8.85 km/h	Up to 8.4km/h	Up to 8.85 km/h	Up to 10 km/h	Up to 9 km/h	Up to 9.65 km/h	Up to 10 km/h	Up to 9.3 km/h
LENGTH & WIDTH ^{3, 11}	L: 119.4 cm W: 52.7 cm	L: 122 cm W: 61.6 cm	L: 115.6 cm W: 57.1 cm	L:115 cm W: 56.5 cm	L: 119.4 cm W: 56.2 cm	L: 119.4 cm W: 59 cm	L: 119.4 cm W: 56.5 cm	3W L: 116.8 cm 3W W: 61 cm 4W L: 118.7 cm 4W W: 61 cm	L: 130.2 cm W: 64.1 cm	3W L: 120.7 cm 3W W: 68.6 cm 4W L: 120.7 cm 4W W: 68.6 cm
SEAT TYPE (STANDARD) ⁵	Compact, Foldable, 43 cm x 43 cm (17" x 17")	Pillow Top Deluxe H-B w/sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Foldable Molded Plastic 46 cm x 43 cm (18" x 17")	Foldable Molded Plastic 46 cm x 43 cm (18" x 17")	CRS Junior Seat 44.25 cm x 43 cm (17.5" x 17")	Foldable Molded Plastic 48 cm x 42 cm (19" x 16.5")	Reclining H-B w/ Sliders & Headrest 46 cm x 46 cm (18" x 18")	Deluxe Reclining H-B w/ Sliders & Headrest 46 cm x 46 cm (18" x 18")	Deluxe Reclining H-B w/ Sliders & Headrest 46 cm x 46 cm (18" x 18")	Adjustable H-B w/Black Vinyl 55.9 cm x 50.8 cm (22" x 20") ²
SEAT-TO-DECK HEIGHT RANGE	38.1 cm - 45.7 cm	47.5 cm - 52.5 cm	37 cm - 42 cm	37 cm - 42 cm	43.8 cm - 47.6 cm	44.5 cm - 48.3 cm	46 cm - 49 cm	41 cm x 44.5 cm	43.2 cm x 53.3 cm	44 cm x 60 cm
SEAT-TO-GROUND HEIGHT RANGE	53.3 cm - 60 cm	58.4 cm - 63.5 cm	51.4 cm - 55.2 cm	52 cm - 55.8 cm	60 cm - 63.8 cm	61 cm - 64 cm	60 cm - 64.8 cm	58.5 cm x 61.5 cm	54 cm - 64 cm	56.6 cm - 72.6 cm
GROUND CLEARANCE ^{3, 12, 14}	8.9 cm at frame 3.2 cm at motor	5 cm at motor	8.8 cm at frame	7.6 cm at frame	8.25 cm at frame	4 cm at rear frame	5.6 cm at motor	8.25 cm at transaxle	10.8 cm at transaxle	5.7 cm at motor
OVERALL WEIGHT (NOT INC. BATTERY) ^{3, 13}	60 kg	82 kg	59.4 kg	62.4 kg	84.8 kg	68.9 kg	81.2 kg	3W: 70.7 kg 4W: 75 kg	112 kg	3W: 86.4 kg 4W: 92 kg
WEIGHT OF HEAVIEST PIECE ³	25.9 kg (front section)	34 kg (rear section)	25.8 kg (front section)	27.7 kg (rear section)	29.5 kg (front section)	50.8 kg (main frame - does not disassemble)	35.2 kg (rear section)	3W: 27.5 kg 4W: 32 kg (front section)	37.6 kg (front section)	3W: 61.7 kg 4W: 67.3 kg (main frame - does not disassemble)
SUSPENSION	Front & Rear	Front & Rear	N/A	N/A	Front & Rear	Front & Rear	Front & Rear	3W: Rear 4W: Front & Rear	Front & Rear	N/A
3W/4W OPTION	4W	4W	4W	4W	4W	4 W	4W	3W & 4W	4W	3W & 4W
FRONT TYRES	2.5" x 9" solid	9" x 2.75" solid	3.6" x 10.4"	3.5" x 10"	3.5" x 10"	3" x 10"	3.5" x 10"	3.5" x 10"	4" x 13"	3.6" x 10" solid
REAR TYRES	2.5" x 9" solid	10.75" x 6" solid	3.6" x 10.75" solid	3.5" x 10"	4" x 10"	3" x 10"	3.5" x 10"	4" x 10"	4" x 10"	3.6" x 10.75" solid
BATTERY REQUIREMENTS ⁶	(2) 31AH	(2) 40AH	(2) 31AH	(2) 31AH or (2) 40AH	(2) 31AH or (2) 40AH	(2) 40AH	(2) 40AH	(2) 31AH or (2) 40AH	(2) 40AH	(2) 50AH
BATTERY CHARGER	Off-Board, 3.5A	Off-Board, 3.5A	Off-Board, 3A	Off-Board, 3A (32AH) Off-Board, 5A (40AH)	Off-Board, 3A (32AH) Off-Board, 5A (40AH)	Off-Board, 5A	Off-Board, 5A	Off-Board, 5A	Off-Board, 5A	Off-board, 5A
BATTERY WEIGHT ¹⁰	11 kg	14 kg	10.5 kg	10.5 kg or 14 kg	10.5 kg or 14 kg	14 kg	14 kg	10.5 kg or 14 kg	14 kg	16.6 kg
COLOUR OPTIONS	Grey Street, True Blue	Candy Apple Red	Red	Viper Blue	Viper Blue	Red	Candy Apple Red	Candy Apple Red	Candy Apple Red, Silver	Garnet Red, Ocean Blue

1 All shroud panels included. **2** Varies with user weight, terrain type, battery amp- hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%. **3** Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. **4** 12 AH battery pack weighs 8.8 kg (19.5 lbs.). **5** Durable solid plastic. **6** Deep cycle (AGM or gel cell type recommended) excluding Lithium-ion models. **7** 18 AH battery pack weighs 13.4 kg (29.5 lbs.). **8** Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. **9** Up to 14 miles per charge on 18 AH batteries. **10** Battery weight may vary based on manufacturer. **11** Width of base (not including seat). **12** Does not account for

anti-tip clearance. **13** Includes standard seat weight. **14** At midframe **15** Range tested with 31-36 AH Batteries. Standard larger battery size would increase range. **16** 6.9 cm x 35.5 cm (2.75" x 14") Rim Only **17** Varies with user weight, terrain type, battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%. **NOTE** The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



END-TO-END | Scooter Specifications Continued...



	PATHRIDER® 130 XL	PATHRIDER® 140 XL	CELEBRITY MFP	OUTBACK®	SPORT RIDER
RANGE PER CHARGE (MAX) ^{2, 8}	Up to 25 km	Up to 35 km	Up to 43.5 km	Up to 35 km	Up to 35 km
MAXIMUM WEIGHT CAPACITY	158.7 kg	181.4 kg	159 kg	158.7 kg	181.4 kg
TURNING RADIUS ³	170 cm	210 cm right 208 cm left	158 cm	243 cm	158 cm
MAXIMUM SPEED ¹⁷	Up to 10 km/h	Up to 10 km/h	Up to 10 km/h	Up to 10 km/h	Up to 10 km/h
LENGTH & WIDTH ^{3, 11}	L: 136 cm W: 63.6 cm	L: 142.8 cm W: 69.2 cm	L: 145 cm W: 70.7 cm	L: 156 cm W: 83 cm	L: 165 cm W: 78 cm
SEAT TYPE (STANDARD) ⁵	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Pillow Top Deluxe H-B w/Sliders 49.5 cm x 44.5 cm (19.5" x 17.5")	Captains High-Back Adjustable Seat 51 cm x 56 cm (20" x 22")
SEAT-TO-DECK HEIGHT RANGE	42 cm - 47 cm	47 cm	45.2 cm	50 cm - 53 cm	47.5 cm
SEAT-TO-GROUND HEIGHT RANGE	56.7 cm	66.7 cm - 70.5 cm	66.7 cm - 70.5 cm	76 cm - 80 cm	24.8 cm
GROUND CLEARANCE ^{3, 12, 14}	9 cm at transaxle	13.3 cm at midframe	8.8 cm at motor	14 cm at midframe	12.7 cm at shroud
OVERALL WEIGHT (NOT INC. BATTERY) ^{3, 13}	99.8 kg	117.7 kg	104.3 kg	124 kg	115 kg
WEIGHT OF HEAVIEST PIECE ³	34 kg (rear section)	93 kg (main frame)	94 kg (main frame - does not disassemble)	99 kg (main frame)	92.3 kg (main frame)
SUSPENSION	Front & Rear	Front & Rear	Front & Rear	Front & Rear	Front & Rear
3W/4W OPTION	4W	4W	4W	4W	3W
FRONT TYRES	4" x 13"	4" x 13"	13"	14.5" x 6"	2.75" x 19" ¹⁶
REAR TYRES	4" x 13"	4" x 13"	32.8 cm (13")	14.5" x 6"	2.75" x 19" ¹⁶
BATTERY REQUIREMENTS ⁶	(2) 55AH	(2) 70AH	(2) 70AH	(2) 100AH	(2) 70AH
BATTERY CHARGER	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A
BATTERY WEIGHT ¹⁰	16.8 kg	23.45 kg	22 kg ⁵	27.9 kg	23.7 kg
COLOUR OPTIONS	Candy Apple Red, Mercury Silver	Candy Apple Red, Yellow, Silver	Red, Blue	Stone Grey	Black



Stealth Products® Cushions feature patented Coolcore® technology, a chemical free, four-way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. Our cushions also use a colour guide that identifies the suggested function family by a unique patch.

FEATURES

- Chemical free
- Regulates temperature
- Manages and distributes your body heat

CUSHIONS.STEALTHPRODUCTS.COM/

Learn about our Stealth Products range of cushions featuring our exclusive Coolcore® Technology.



Scan here to access the
Stealth Cushions Catalogue



30+ UPF





Dry



Chemical
Free



Breathable
Comfort



***Refer to previous pages for footnotes

Jazzy® POWER CHAIRS



Go Chair

Jazzy® PASSPORT



Embrace the ease of lightweight travel with the Jazzy® Passport! The compact design effortlessly folds in a few simple steps, making it a breeze to transport. With convenient storage capabilities, the Jazzy® Passport ensures hassle-free traveling.

 Range per Charge: <15.4 km

 Weight Capacity: 113.3 kg

 Turning Radius: 61.5 cm

GOCHAIR®



The Go Chair® NG is #1 for portability! With Pride's feather-touch disassembly, and a tight turning radius, the Go Chair® NG is ideal for indoor use and smooth surfaces, allowing you to enjoy greater independence.

 Range per Charge: <14 km

 Weight Capacity: 136 kg

 Turning Radius: 64.5 cm

Jazzy® SELECT 6



Experience the sporty and sophisticated design of the Jazzy® Select 6 power wheelchair. It provides powerful performance and stability with its superior suspension system. Featuring an adjustable high-back seat with manual recline.

Range per Charge: <24.1 km

Weight Capacity: 136 kg

Turning Radius: 55.5 cm

Jazzy® EVO 613Li

NEW



The Jazzy® EVO 613 provides the freedom to navigate small spaces and tight corners indoors with a narrow base width of only 56cm. The Jazzy® EVO 613Li uses a lithium-ion battery that provides longer life and more convenient charging.


Range per Charge: <27.5 km

Weight Capacity: 136 kg

Turning Radius: 55.8 cm

DOWNLOAD OUR PQ-365 APP
FOR EVERYTHING PRIDE AND QUANTUM

PQ-365



Download on the
App Store

GET IT ON
Google Play

Jazzy® air₂



Take your life to new heights with the Jazzy Air® 2! With a simple switch, this power wheelchair elevates 30.4 cm in just 11 seconds. Allowing you to safely drive at speeds of up to 6.4 km/h while elevated, enabling face-to-face social interactions.

Range per Charge: <27 km

Weight Capacity: 136 kg

Turning Radius: 61.5 cm

Jazzy® 623 2.0



The Jazzy® 623 2.0 is engineered to meet the performance needs of the most active user. Equipped with an adjustable high-back seat and adjustable armrests, as well as powerful motors, this power wheelchair ensures exceptional performance.

 Range per Charge: <24.9 km

 Weight Capacity: 136 kg

 Turning Radius: 50.8 cm



Jazzy Air 2

Jazzy® ELITE HD



Experience the ultimate heavy-duty electric wheelchair with the Jazzy® Elite HD. This power wheelchair excels on various terrains, providing superior performance and control with its dual-motor, front-wheel drive technology, that offers a tight turning radius.

 Range per Charge: <15.7 km

 Weight Capacity: 204.1 kg

 Turning Radius: 61.9 cm

Jazzy® 1450



With a robust and heavy-duty design, the Jazzy® 1450 ensures powerful performance with front-wheel drive and powerful motors. This bariatric chair offers a 24" - 32" seat size, allowing you to live life on your terms.

 Range per Charge: <25.2 km

 Weight Capacity: 273.1 kg

 Turning Radius: 67.3 cm

END-TO-END | Powerchair Specifications



	JAZZY® PASSPORT	GO CHAIR® NG	JAZZY® SELECT 6	JAZZY® EVO 613 LI	JAZZY® AIR 2	JAZZY® 623 2.0	JAZZY® ELITE HD	JAZZY® 1450
RANGE PER CHARGE (MAX) ^{9,5}	Up to 15.4 km	Up to 14km	Up to 24.1 km	Up to 27.5 km	Up to 27 km	Up to 24.9 km	Up to 15.7 km	Up to 25.2 km
MAXIMUM WEIGHT CAPACITY	113.3 kg	136 kg	136 kg	136 kg	136 kg	136 kg	204.1 kg ⁶	273.1 kg
TURNING RADIUS ¹	61.5 cm	64.5 cm	55.5 cm	55.8 cm	61.5 cm	50.8 cm	61.9 cm	67.3 cm
MAXIMUM SPEED ³	Up to 5.7 km/h	Up to 5.9 km/h	Up to 6.4 km/h	Up to 7.8 km/h	Up to 6.4 km/h	Up to 9.6 km/h	Up to 6.4 km/h	Up to 7.2 km/h
LENGTH & WIDTH ¹	L: 101.6 cm W: 59.6 cm	L: 83.8 cm W: 57 cm	L: 90 cm W: 59.6 cm	L: 104 cm W: 55.8 cm	L: 106 cm W: 61.5 cm	L: 90 cm W: 60.9 cm	L: 90 cm W: 61.5 cm	L: 97 cm W: 74.2 cm
WEIGHT OF HEAVIEST PIECE ⁴	24.4 kg	16.3 kg	44.1 kg	42.6 kg	91.6 kg ⁸	62.3 kg	44.9 kg	88.4 kg
GROUND CLEARANCE ¹	7.6 cm at motor	4.5 cm at frame	4.4 cm at motor mount	7.6 cm at motors	3.6 cm at motor mount	7.6 cm at frame	7.6 cm at footboard mount 3.8 cm at anti-tip wheel	7.6 cm at footboard mount 3.8 cm at anti-tip wheel
SEAT SIZE (STANDARD)	41.9 cm x 44.4 cm (16.5"x 17.5")	43.1 cm x 45.7 cm (17"x 18")	45.7 cm x 45.7 cm (18"x 18")	40.6 cm x 40.6 cm (16"x16") 45.7 cm x 45.7 cm (18"x18") 50.8 cm x 45.7 cm (20"x18")	40.6 cm x 45.7 cm (16"x 18") 45.7 cm x 45.7 cm (18"x 18") 50.8 cm x 45.7 cm (20"x 18")	45.7 cm x 45.7 cm (18"x 18") 50.8 cm x 45.7 cm (20"x 18")	45.7 cm x 45.7 cm (18"x 18") 50.8 cm x 50.8 cm (20"x 20") 55.8 cm x 50.8 cm (22"x 20")	60.9 cm to 81.2 cm (24" to 32")
SEAT WEIGHT (STANDARD)	N/A (seat is integrated with the base) (sling)	13.6 kg	19.9 kg (comfort, high-back)	22.7 kg	23 kg (high-back)	19.9 kg (high-back)	23.5 kg (high-back)	25.8 kg (high-back)
HIGH-BACK SEAT-TO-FLOOR RANGE	N/A	N/A	53.9 cm - 59 cm	50.2 cm x 54 cm	53.9 cm - 84 cm	52 cm - 59.6 cm	50.8 cm - 55.8 cm	52 cm - 57 cm
SPECIALITY SEATING OPTIONS	N/A	N/A	N/A	N/A	N/A	Synergy® Seating ^A Tru-Balance ^B	N/A	Synergy® Seating ^A
ELECTRONICS (STANDARD)	50A, PG VSI	40A, PG GC3	40A, PG GC3	40A, PG GC3	40A, VR2	90A, PG VR2	60A, Dynamic Linx	90A, PG VR2
SPECIALTY CONTROLS	N/A	N/A	N/A	N/A	N/A	Yes	N/A	Yes
DRIVE TRAIN	2x 180W brushed motors Rear Drive	Two-motor, rear-wheel drive	Two-motor, in-line, mid-wheel drive	Two motor, in-line, mid-wheel drive	Two-motor, in-line, mid-wheel drive	Two-motor, in-line, mid-wheel drive ⁶	Two-motor, front-wheel drive	Two-motor, front-wheel drive
BRAKING SYSTEM	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical	Regenerative and Electro-mechanical
SUSPENSION TYPE	Limited	Limited	Active-Trac® (ATX)	Active-Trac® (ATX)	Active-Trac® with Mid-Wheel 6 technology	Active-Trac® (ATX)	Limited	Limited
FRONT WHEELS	8" solid castors	6" solid castors	6" solid castors	5" solid	6" solid castors	6" solid castors	3" solid anti-tips	4" solid anti-tips
DRIVE WHEELS	12" pneumatic	9" solid	10" solid	13" solid	10" solid	14" pneumatic	14" solid	14" solid
REAR WHEELS	N/A	3" solid anti-tips	6" solid castors	5" casters	6" solid castors	8" solid castors	6" solid castors	9" castors
BATTERY REQUIREMENTS	24 V, 18 AH Lithium-ion (1)	12 V, 18 AH (2)	12 V, 35 AH (2)	(1) 25V Lithium-iron phosphate	12 V, 40 AH (2)	12 V, 55 AH (2)	12 V, 55 AH (2)	12 V, 70 AH (2)
BATTERY CHARGER	2A, off-board, Lithium	2A, off-board	3.5A, off-board	3.5A, Lithium Charger	5A, off-board	8A, off-board	5A, off-board	8A, off-board
BATTERY WEIGHT ⁷	3.9 kg	14 kg battery box	10.5 kg each	8.4 kg	14 kg each	16.8 kg each	16.8 kg each	23.4 kg each
COLOUR OPTIONS	Pewter	Ruby Red, Pearl White	Ruby Red	Black, Robins Egg Blue, Cardinal Red	Black Oynx, Pearl White, Pink Topaz, Garnet Red	Candy Apple Red	Jazzy Red, Jazzy Blue	Candy Apple Red, Viper Blue

¹ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. ² Without front riggings. ³ Varies with user weight, terrain, battery condition and tire pressure. ⁴ Does not include seat or battery weights. ⁵ Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, - 5%. ⁶ Maximum weight capacity is reduced to 400 lbs. (181 kg) on 46 cm (18") deep seat configurations. ⁷ Battery weight may vary based on manufacturer. ⁸ Seat weight included (because tools are required to remove the seat). Battery

weight not included. ⁹ Varies with user weight, terrain, battery type, battery condition and tire pressure. ^{*}International configuration shown **SEATING: A** For Synergy® Seating Contact your nearest Pride® stocklist for more information **B** Quantum® dealers only **NOTE** The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



PRIDE RAMP

SINGLE FOLD




Single Fold Ramps offer a reliable solution for transitioning curbs, thresholds, and other obstacles. They are 75 cm wide and are available in several different lengths. These high-traction ramps are portable and lightweight, making them ideal for home and travel.



Weight Capacity: 272 kg

RUBBER THRESHOLD



Rubber Threshold Ramps provide a safe, durable means of navigating uneven household surfaces, including small steps and doorways. A hidden dovetail locking system allows for easy installation and incremental ramp height adjustment of up to 10.16cm (4"). A modular design, slip-resistant and maintenance free material makes our Rubber Threshold Ramps a practical and affordable home accessibility solution.




Weight Capacity: 385 kg

ACCESSORIES



Sun Canopy & Flag



Rear View Mirror



Rear Pride Scooter Bag



Saddle Pride Bag



Walking Frame Holder



Oxygen Tank Holder



Rear Pod



XLR USB Charger



Mobile Phone Holder



Foldable Cup Holder



Tel: (03) 8770 9600

Visit pridemobility.com.au for further details



VIC Head Office: 20-24 Apollo Dr. Hallam, VIC 3803

NSW Distribution Centre: 19 Newton Rd. Wetherill Park, NSW 2164

QLD Distribution Centre: Unit 3, 29 McCotter St. Acacia Ridge, QLD 4110

NZ Distribution Centre: Factory 5/208 Swanson Rd. Henderson, Auckland NZ 0610

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.