

The most advanced Jazzy® ever

The advanced comfort and versatility of the Jazzy® EVO 613 Li provides the freedom to navigate small spaces and tight corners. Its narrow base width of only 56cm makes it an ideal powerchair for moving through your home. The Jazzy® EVO 613 Li uses a lightweight lithium-ion battery that provides longer life and 25% faster charging. The high-back swivel seat offers convenience with an additional option of a swing away controller mount.











Optional Swing Away Controller Mount



Mid-Wheel 6® Technology

RANGE PER CHARGE (MAX) ^{1,3}	Up to 27.5 km
MAXIMUM WEIGHT CAPACITY	136 kg
TURNING RADIUS ²	55.8 cm
MAXIMUM SPEED ¹	7.8 km/h
LENGTH & WIDTH ²	L: 104 cm W: 55.8 cm
WEIGHT OF HEAVIEST PIECE ⁴	42.6 kg
GROUND CLEARANCE ²	7.6 cm at motors
SEAT SIZE	L: 45.7 cm (18 - 20") W: 45.7 cm (18)
STANDARD 18" X 18" SEAT WEIGHT	22.7 kg
HIGH-BACK SEAT-TO-FLOOR RANGE	50.2 cm - 54 cm
SPECIALTY SEATING OPTIONS	N/A

- Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
- 2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
- 4. Does not indluce seat or battery weights.
- 5. Battery weight may very based on manufacturer.

ELECTRONICS	40A, PG GC3
SPECIALTY CONTROLS	N/A
DRIVE TRAIN	Two motor, in-line, mid-wheel drive
BRAKING SYSTEM	Regenerative and electro-mechanical
SUSPENSION TYPE	Active-Trac® (ATX)
FRONT WHEELS	5" solid castors
DRIVE WHEELS	13" solid
REAR WHEELS	5" solid castors
BATTERY REQUIREMENTS	25.6 V, 32.4 AH Lithium-ion (1)
BATTERY CHARGER	3.5A, off-board, Lithium
BATTERY WEIGHT ⁵	8.4 kg
COLOUR OPTIONS	Black, Robins Egg Blue, Cardinal Red

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.



