



|                                       | JAZZY <sup>®</sup> PASSPORT                    | GO CHAIR <sup>®</sup> NG               | JAZZY <sup>®</sup> SELECT 6            | JAZZY <sup>®</sup> EVO 613 LI   | JAZZY <sup>®</sup> AIR 2   | JAZZY <sup>®</sup> 623 2.0  | JAZZY <sup>®</sup> ELITE HD  | JAZZY <sup>®</sup> 1450                               |
|---------------------------------------|--|--|--|---|--|---|--|---|
| RANGE PER CHARGE (MAX) <sup>9,5</sup> | Up to 15.4 km                                  | Up to 14km                             | Up to 24.1 km                          | Up to 27.5 km   | Up to 27 km  | Up to 24.9 km   | Up to 15.7 km  | Up to 25.2 km   |
| MAXIMUM WEIGHT CAPACITY               | 113.3 kg                                       | 136 kg                                 | 136 kg                                 | 136 kg  | 136 kg   | 136 kg  | 204.1 kg <sup>6</sup>  | 273.1 kg  |
| TURNING RADIUS <sup>1</sup>           | 61.5 cm  | 64.5 cm                                | 55.5 cm                                | 55.8 cm   | 61.5 cm  | 50.8 cm   | 61.9 cm  | 67.3 cm   |
| MAXIMUM SPEED <sup>3</sup>            | Up to 5.7 km/h                                 | Up to 5.9 km/h                         | Up to 6.4 km/h                         | Up to 7.8 km/h  | Up to 6.4 km/h   | Up to 9.6 km/h  | Up to 6.4 km/h   | Up to 7.2 km/h  |
| LENGTH & WIDTH <sup>1</sup>           | L: 101.6 cm W: 59.6 cm                         | L: 83.8 cm W: 57 cm                    | L: 90 cm W: 59.6 cm                    | L: 104 cm W: 55.8 cm  | L: 106 cm W: 61.5 cm   | L: 90 cm W: 60.9 cm   | L: 90 cm W: 61.5 cm  | L: 97 cm W: 74.2 cm                                   |
| WEIGHT OF HEAVIEST PIECE <sup>4</sup> | 24.4 kg  | 16.3 kg                                | 44.1 kg                                | 42.6 kg   | 91.6 kg <sup>8</sup>   | 62.3 kg   | 44.9 kg  | 88.4 kg   |
| GROUND CLEARANCE <sup>1</sup>         | 7.6 cm at motor                                | 4.5 cm at frame                        | 4.4 cm at motor mount                  | 7.6 cm at motors  | 3.6 cm at motor mount  | 7.6 cm at frame   | 7.6 cm at footboard mount<br>3.8 cm at anti-tip wheel  | 7.6 cm at footboard mount<br>3.8 cm at anti-tip wheel |
| SEAT SIZE                             | 41.9 cm x 44.4 cm<br>(16.5"x 17.5")            | 43.1 cm x 45.7 cm<br>(17"x 18")        | 45.7 cm x 45.7 cm<br>(18"x 18")        | 40.6 cm x 40.6 cm (16"x16")<br>45.7 cm x 45.7 cm (18"x18")<br>50.8 cm x 45.7 cm (20"x18") | 40.6 cm x 45.7 cm (16"x 18")<br>45.7 cm x 45.7 cm (18"x 18")<br>50.8 cm x 45.7 cm (20"x 18") | 45.7 cm x 45.7 cm<br>(18"x 18")<br>50.8 cm x 45.7 cm<br>(20"x 18")    | 45.7 cm x 45.7 cm (18"x 18")<br>50.8 cm x 50.8 cm (20"x 20")<br>55.8 cm x 50.8 cm (22"x 20") | 60.9 cm to 81.2 cm<br>(24" to 32")                    |
| 18' X 18" STANDARD SEAT WEIGHT        | N/A (seat is integrated with the base) (sling) | 13.6 kg                                | 19.9 kg (comfort, high-back)           | 22.7 kg   | 23 kg (high-back)  | 19.9 kg (high-back)   | 23.5 kg (high-back)  | 25.8 kg (high-back)                                   |
| HIGH-BACK SEAT-TO-FLOOR RANGE         | N/A  | N/A                                    | 53.9 cm - 59 cm                        | 50.2 cm - 54 cm   | 53.9 cm - 84 cm  | 52 cm - 59.6 cm   | 50.8 cm - 55.8 cm  | 52 cm - 57 cm   |
| SPECIALITY SEATING OPTIONS            | N/A  | N/A                                    | N/A                                    | N/A   | N/A  | Synergy <sup>®</sup> Seating <sup>A</sup><br>Tru-Balance <sup>B</sup> | N/A  | Synergy <sup>®</sup> Seating <sup>A</sup>             |
| STANDARD ELECTRONICS                  | 50A, PG VSI                                    | 40A, PG GC3                            | 40A, PG GC3                            | 40A, PG GC3   | 40A, VR2   | 90A, PG VR2   | 60A, Dynamic Linx  | 90A, PG VR2   |
| SPECIALTY CONTROLS                    | N/A  | N/A                                    | N/A                                    | N/A   | N/A  | Yes   | N/A  | Yes   |
| DRIVE TRAIN                           | 2x 180W brushed motors<br>Rear Drive           | Two-motor,<br>rear-wheel drive         | Two-motor, in-line,<br>mid-wheel drive | Two motor, in-line,<br>mid-wheel drive  | Two-motor, in-line,<br>mid-wheel drive   | Two-motor, in-line,<br>mid-wheel 6 <sup>®</sup>                       | Two-motor,<br>front-wheel drive  | Two-motor,<br>front-wheel drive                       |
| BRAKING SYSTEM                        | Regenerative and<br>Electro-mechanical         | Regenerative and<br>Electro-mechanical | Regenerative and<br>Electro-mechanical | Regenerative and<br>Electro-mechanical  | Regenerative and<br>Electro-mechanical   | Regenerative and<br>Electro-mechanical                                | Regenerative and<br>Electro-mechanical   | Regenerative and<br>Electro-mechanical                |
| SUSPENSION TYPE                       | Limited  | Limited                                | Active-Trac <sup>®</sup> (ATX)         | Active-Trac <sup>®</sup> (ATX)  | Active-Trac <sup>®</sup> with<br>Mid-Wheel 6 technology                                      | Active-Trac <sup>®</sup> (ATX)  | Limited  | Limited   |
| FRONT WHEELS                          | 8" solid castors                               | 6" solid castors                       | 6" solid castors                       | 5" solid  | 6" solid castors   | 6" solid castors  | 3" solid anti-tips   | 4" solid anti-tips                                    |
| DRIVE WHEELS                          | 12" pneumatic                                  | 9" solid                               | 10" solid                              | 13" solid   | 10" solid  | 14" pneumatic   | 14" solid  | 14" solid   |
| REAR WHEELS                           | N/A  | 3" solid anti-tips                     | 6" solid castors                       | 5" casters  | 6" solid castors   | 8" solid castors  | 6" solid castors   | 9" castors  |
| BATTERY REQUIREMENTS                  | 24 V, 18 AH Lithium-ion (1)                    | 12 V, 18 AH (2)                        | 12 V, 35 AH (2)                        | 25.6 V, 32.4 AH Lithium-ion (1)   | 12 V, 40 AH (2)  | 12 V, 55 AH (2)   | 12 V, 55 AH (2)  | 12 V, 70 AH (2)                                       |
| BATTERY CHARGER                       | 2A, off-board, Lithium                         | 2A, off-board                          | 3.5A, off-board                        | 3.5A, off-board, Lithium  | 5A, off-board  | 8A, off-board   | 5A, off-board  | 8A, off-board   |
| BATTERY WEIGHT <sup>7</sup>           | 3.9 kg   | 14 kg battery box                      | 10.5 kg each                           | 8.4 kg  | 14 kg each   | 16.8 kg each  | 16.8 kg each   | 23.4 kg each  |
| COLOUR OPTIONS                        | Pewter   | Ruby Red, Pearl White                  | Ruby Red                               | Black, Robins Egg Blue,<br>Cardinal Red   | Black Oynx, Pearl White,<br>Pink Topaz, Garnet Red   | Candy Apple Red   | Jazzy Red, Jazzy Blue  | Candy Apple Red, Viper Blue                           |

<sup>1</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. <sup>2</sup> Without front riggings. <sup>3</sup> Varies with user weight, terrain, battery condition and tire pressure. <sup>4</sup> Does not include seat or battery weights. <sup>5</sup> Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, - 5%. <sup>6</sup> Maximum weight capacity is reduced to 400 lbs. (181 kg) on 46 cm (18") deep seat configurations. <sup>7</sup> Battery weight may vary based on manufacturer. <sup>8</sup> Seat weight included (because tools are required to remove the seat). Battery weight not included. <sup>9</sup> Varies with user weight, terrain, battery type, battery condition and tire pressure. \*International configuration shown **SEATING: A** For Synergy<sup>®</sup> Seating Contact your nearest Pride<sup>®</sup> stocklist for more information **B** Quantum<sup>®</sup> dealers only **NOTE** The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.